

Colorimetric PBG Assay (Ehrlich)

COLORIMETRIC PORPHOBILINOGEN ASSAY (EHRlich)

Modified Ehrlich Reagent

1. Dissolve 1.0 gm p-dimethyl aminobenzaldehyde in 30mL glacial acetic acid.
2. Add 8.0 mL 70% perchloric acid.
3. Bring to 50 mL with glacial acetic acid.
4. Store in the dark at 4°C for about 2 weeks.

Porphobilinogen (PBG) Assay

1. Mix 2.0mL sample with 2.0mL fresh modified Ehrlich reagent.
2. Measure absorbance at 555nm.

Standard PBG

1. Prepare 0.54mg/mL (2.2mM) PBG stock in 0.1M Tris pH 7.65.
2. (Store at -20°C until needed.)
3. Dilute 1:1000 in H₂O.
4. Mix 2mL with 2mL fresh modified Erlich reagent.
5. Measure absorbance at 555nm.
 - 1mM PBG = 61 abs units at A555.

Revision #2

Created 2023-12-01 21:55:39 UTC by Elliot Francis

Updated 2023-12-04 16:29:06 UTC by Elliot Francis