

FECH Activity Assay

FERROCHELATASE (FECH)

Sample preparation

1. Suspend ~50- μ L mammalian cell pellet in 400 μ L TGD buffer (Tris buffered glycerol with dithiothreitol (DTT), made by dissolving 2mL or 2.52 g glycerol and 1.5mg DTT in 8.0 mL 20mM Tris pH 8.0)
2. Sonicate (homogenize) while in an ice bath at the lowest practicable power setting for 3 cycles x 5 sec at 50% duty (2.5 sec on, 2.5 sec off)
3. Determine the protein content and dilute to 1 μ g protein/ μ L with more TGD

Ferrochelatase reaction

1. Prepare three 50- μ L aliquots of the cell preparation, two live and one inactivated in boiling water for 10 minutes (as control for non-enzymatic product formation).
2. Prepare the incubation buffer containing 160mM Tris pH 8.0, 40mM Bicine pH 8.0, 10mg/ml Tween20 and 0.38mg/mL palmitic acid),
3. Mix 150 μ L incubation buffer and 25 μ L zinc substrate (1mM aqueous Zn acetate).
4. Mix each 50- μ L aliquot of cell preparation with the 150 μ L incubation buffer plus zinc substrate and pre-incubate for 5 minutes at 37°C.
5. Then add 25 μ L of mesoporphyrin IX substrate (250 μ M in 160mM Tris pH 8.0, 40mM Bicine pH8.0, 2mg/ml Tween20).
6. Incubate the mixture for 30 min at 37°C.
7. Add 750 μ L stop reagent (270 μ M ethylenediaminetetraacetic acid in a mixture containing dimethylsulfoxide-methanol, 30/70 by volume respectively).
8. Cool on ice for 15-20 min.
9. Centrifuged at 1500xg for 10 min at room temperature.

Quantitation of product Zn mesoporphyrin IX

1. Inject 10 μ L of supernatant solution of porphyrins into a Waters Acquity UPLC (ultra performance liquid chromatography) system, which includes a binary solvent manager, sample manager, fluorescence detector (FLR), column heater and an Acquity UPLC BEH C18, 1.7 μ M, 2.1 x 100 mm column.
2. Set the FLR for zinc mesoporphyrin IX (ZnMeso) at 406 nm excitation and 578 nm emission.
3. Quantify the ferrochelatase product relative to a standard ZnMeso solution, also in the stop reagent.
4. Keep the sample chamber dark and at ambient temperature.
5. Solvent A is 0.2% aqueous formic acid while Solvent B was 0.2% formic acid in methanol.

6. Set the flow rate at 0.40 mL per minute at 60°C for the total run time of 7 min.
 7. Use the following successive gradient settings for run time in minutes versus A:
 - 0.0, 80%
 - 2.5, 1%
 - 4.5, 1%
 - 5.0, 80%
 8. Keep all solvent gradients are linear.
-

Revision #2

Created 2023-12-01 21:55:16 UTC by Elliot Francis

Updated 2023-12-04 16:26:14 UTC by Elliot Francis