

Primer Sequences

These are the sequences of the primers that we have aliquots of here in the core for our users to make use of.

M13F(-21)	GTAAAACGACGGCCAGT 3'
M13R	CAGGAAACAGCTATGAC 3'
T7prom	TAATACGACTCACTATAGG 3'
T7term	GCTAGTTATTGCTCAGCGG 3'
T3prom	ATTAACCCTCACTAAAGGG3'
SP6	GATTTAGGTGACACTATAG 3'
CMV-F	CGCAAATGGGCGGTAGGCGTG 3'
CMV-R	AGTAGGAAAGTCCCGTAAGG
EGFP-C	CATGGTCCTGCTGGAGTTCGTG
EGFP-N	CGTCGCCGTCCAGCTCGACCAG

All users are welcome to come get an aliquot of 50uL of 2uM primer dilutions for use on samples here through our core. Please do NOT re-freeze the aliquots. They will store fine in the fridge. Repeated freeze/thaw cycles of unconcentrated primers can be problematic is Sanger sequencing.

Revision #1

Created 2024-03-18 17:53:56 UTC by Derek Warner

Updated 2024-03-18 17:59:35 UTC by Derek Warner